


Status	Released	ID/Version/Description	MetaT002/00000001/Meta Basic MBR
Valid from		Standard batch size	1 (0.5 - 1.5) kg
		Estimated yield	1 kg

Approval		
Status	Text	
Master Batch Record 10#MetaT002#Meta Basic MBR#Meta0001#---##00000001		
Draft	Object was created	Date - User: 04-09-2024 13:18:47 CEST - 110 / MBR_Designers
In circulation	Put MBR in circulation	Date - User: 05-09-2024 15:20:38 CEST - 110 / MBR_Designers
Comment: In circulation		
In circulation	MBR approval by technical support team	Date - User: 05-09-2024 15:52:00 CEST - 180 / PASX_Engineering
Comment: Approved		
Draft	MBR rejection by technical support team	Date - User: 05-09-2024 15:53:14 CEST - 180 / PASX_Engineering
Comment: Rejected because Valerie said it was horrible!!!!!!!!		
In circulation	Put MBR in circulation	Date - User: 05-09-2024 15:55:27 CEST - 110 / MBR_Designers
Comment: In circ again		
In circulation	MBR approval by technical support team	Date - User: 09-09-2024 11:30:10 CEST - 180 / PASX_Engineering
Released	MBR approval by technical support team	Date - User: 09-09-2024 16:02:54 CEST - 227 / Valerie.Bogaerts

Basic operations	
Basic operation	
BO1 - Weighing	Active: yes
BO2 - Tableting	Active: yes

Summing bill of input material				
Material no.	Material name	Quantity	Use type	Number of positions
MetaCS01	Cornstarch	0.88 kg	Raw material	1
MetaGL01	Glucose	0.11 kg	Compensating material	1
MetaSA01	Meta Salycilic acid	0.01 kg	Active ingredient with potency and compensating material	1

Summing bill of target material				
Material no.	Material name	Quantity	Use type	Number of positions
Meta0001	Meta Aspirin	1 kg	Raw material	1



Version history		
Version	Status	Text
Master Batch Record 10#MetaT002#Meta Basic MBR#Meta0001#---##00000001		
00000001	Draft	Object was created Date - User: 04-09-2024 13:18:47 CEST - 110 / MBR_Designers
00000001	In circulation	Put MBR in circulation Comment: In circulation Date - User: 05-09-2024 15:20:38 CEST - 110 / MBR_Designers
00000001	In circulation	MBR approval by technical support team Comment: Approved Date - User: 05-09-2024 15:52:00 CEST - 180 / PASX_Engineering
00000001	Draft	MBR rejection by technical support team Comment: Rejected because Valerie said it was horrible!!!!!!! Date - User: 05-09-2024 15:53:14 CEST - 180 / PASX_Engineering
00000001	In circulation	Put MBR in circulation Comment: In circ again Date - User: 05-09-2024 15:55:27 CEST - 110 / MBR_Designers
00000001	In circulation	MBR approval by technical support team Date - User: 09-09-2024 11:30:10 CEST - 180 / PASX_Engineering
00000001	Released	MBR approval by technical support team Date - User: 09-09-2024 16:02:54 CEST - 227 / Valerie.Bogaerts

Active: yes **Reconciliation:** no
Start condition: Predecessors started


Bill of input material											
Item	ERP pos.	Prio.	Material no.	Material name	Optional	Use type	Quantity	Scrap (%)	Quantification type	Manufacturing date	Target material
0010	0	1	MetaSA01 Active: yes	Meta Salycilic acid	no	Active ingredient with potency and compensating material	0.01 kg	0	Weighing	First take-out	PB_SA01
0020	0	2	MetaGL01 Active: yes	Glucose	no	Compensating material	0.11 kg	0	Weighing	No	PB_GL01
0030	0	3	MetaCS01 Active: yes	Cornstarch	no	Raw material	0.88 kg	0	Weighing	No	PB_CS01

Bill of target material											
Item	Material no.			Material name	Optional	Quantity	Usage sign	Material flow	Hand over	Quantification type	Target BO
0040	PB_SA01 Active: yes Container groups			Polybag MetaSA01	no	0.01 kg		Order	no	Manual	BO2
	Number of containers		Target quantity per container		Default variant		Container group ID				
	1				yes		MetaPoly				
0050	PB_GL01 Active: yes Container groups			Polybag MetaGL01	no	0.11 kg		Order	no	Manual	BO2
	Number of containers		Target quantity per container		Default variant		Container group ID				
	1				yes		MetaPoly				
0060	PB_CS01 Active: yes Container groups			Polybag MetaSC01	no	0.88 kg		Order	no	Manual	BO2
	Number of containers		Target quantity per container		Default variant		Container group ID				
	1				yes		MetaPoly				

Production unit		
ID	Description	Production area
PRODUCTION	Production	Packaging Pharmacy

Bill of equipment			
Pos.	Equipment class	Type	Equipment
1	Workroom	MetaWB / Meta Weighing Booth	--- / --- Active: yes

Active: yes **Reconciliation:** no
Start condition: Predecessors started



Detailed BO				
Level		Step		
1		(1) MetaT002/BO1/CHKWB Allocate Weighing Booth Active: yes Optional: no (EQM Allocation) Subtype: None		
			Actual value	Date
				Signature
		Check and allocation of the weighing booth		
		Attributes		
		Identification type: Scan	Allowed status formula:	Reuse equipment: no
		Equipment item: 1: Workroom MetaWB Meta Weighing Booth		
1		(1) MetaT002/BO1/GMPCHK GMP Check Active: yes Optional: no (Common BF) Subtype:		
		Check if material package is accurate		
	1.0	SOP - Read SOP Active: yes Subtype: None	Training Doc (PAS-X v3.1.8 Basic Training.pdf)	---
		Set value	Actual value	Entry / Signature
		Training Doc (PAS-X v3.1.8 Basic Training.pdf)		
	1.1	MATCHK - Material Check Active: yes Subtype: None	Material OK	CQA
		10 Check the material for apparent damages		
		Set value	Actual value	Entry / Signature
		Material OK		
		Date / signature for basic function		
		User identification:	No identification	
		User ident. group:	-----	
		Verification signature:	no	

Active: yes Reconciliation: no
Start condition: Predecessors started

1	(1) MetaT002/BO1/TO1 Take-out Active: yes Optional: no (Take-out) Subtype: None								
	<table><tr><td></td><td>Actual value</td><td>Date</td><td>Signature</td></tr><tr><td></td><td></td><td></td><td></td></tr></table>		Actual value	Date	Signature				
	Actual value	Date	Signature						
	Material take-out and dispensing								
	Attributes								
	All input materials: yes Ignore batch status: no REM unpacking: no								
	Execution mode: Manual								

1	(1) MetaT002/BO1/DA Release Weighing Booth Active: yes Optional: no (EQM Deallocation) Subtype: None								
	<table><tr><td></td><td>Actual value</td><td>Date</td><td>Signature</td></tr><tr><td></td><td></td><td></td><td></td></tr></table>		Actual value	Date	Signature				
	Actual value	Date	Signature						
	Release Weighing Booth								
	Attributes								
	Identification type: Scan Allowed status formula: Reuse equipment: no								
	Equipment item: 1: Workroom MetaWB Meta Weighing Booth								

Active: yes Reconciliation: no
Start condition: Predecessors started


Bill of input material											
Item	ERP pos.	Prio.	Material no.	Material name	Optional	Use type	Quantity	Scrap (%)	Quantification type	Manufacturing date	Target material
0070	0	1	PB_SA01 Active: yes	Polybag MetaSA01	no	Raw material	0.01 kg	0	Take over	No	
0080	0	1	PB_GL01 Active: yes	Polybag MetaGL01	no	Raw material	0.11 kg	0	Take over	No	
0090	0	1	PB_CS01 Active: yes	Polybag MetaSC01	no	Raw material	0.88 kg	0	Take over	No	

Bill of target material								
Item	Material no.	Material name	Optional	Quantity	Usage sign	Material flow	Hand over	Quantification typeTarget BO
0100	Meta0001	Meta Aspirin	no	1 kg		Stock	no	Weighing---
	Active: yes							
	Container groups							
	Number of containers	Target quantity per container	Default variant		Container group ID			
	1		yes		MetaConS			


Production unit		
ID	Description	Production area
PRODUCTION	Production	Packaging Pharmacy

Bill of equipment				
Pos.	Equipment class	Type	Equipment	Number of objects
1	Workroom	MetaTab / Meta Tableting	--- / --- Active: yes	1
2	Setup part	MetaSP / Meta set-up parts	--- / --- Active: yes	1

Active: yes Reconciliation: no
Start condition: Predecessors started



Active: yes Reconciliation: no
Start condition: Predecessors started



Detailed BO				
Level Step				
1		(1) MetaT002/BO2/TO2 Material Take Over Active: yes Optional: no (Take-out) Subtype: None		!
			Actual value	Date
				Signature
		Take over materials		
		Attributes		
		All input materials: yes	Ignore batch status: no	REM unpacking: no
		Execution mode: Manual		
1		(1) MetaT002/BO2/AS1 EQM Assembly/Setup Active: yes Optional: no (EQM Assembly/Setup) Subtype: None		
			Actual value	Date
				Signature
		Install Setup Part in tablet press		
		Attributes		
		Identification type: Scan	Allowed status formula:	Implicit allocation: yes
		Status check: Diagram based	Reuse equipment: no	
		Equipment item: 1: Workroom MetaTab Meta Tableting		
		Bill of equipment		
		Pos.	Equipment	Number of objects
		2 Setup part MetaSP / Meta set-up parts	--- / ---	1
		Active: yes		

Active: yes **Reconciliation:** no
Start condition: Predecessors started

1	(1) MetaT002/BO2/IPC In process control Active: yes Optional: no (Common BF)						
	Subtype: Check humidity of the material						
1.0	HUM01 - Check Humidity Active: yes Subtype: None 45.0 (≥ 10.0 ; ≤ 90.0) CPP						
	10 Check humidity with hygrometer device H6643G and enter value						
	<table><tr><th>Set value</th><th>Actual value</th><th>Entry / Signature</th></tr><tr><td>45.0 (≥ 10.0 ; ≤ 90.0)</td><td></td><td></td></tr></table>	Set value	Actual value	Entry / Signature	45.0 (≥ 10.0 ; ≤ 90.0)		
Set value	Actual value	Entry / Signature					
45.0 (≥ 10.0 ; ≤ 90.0)							
	Date / signature for basic function User identification: No identification User ident. group: ----- Verification signature: no						

1	(1) MetaT002/BO2/DEC01 Decision Optional: NO (DECISION)												
	Subtype: NONE												
	<table><tr><th>Actual value</th><th>Date</th><th>Signature</th></tr><tr><td></td><td></td><td></td></tr></table>	Actual value	Date	Signature									
Actual value	Date	Signature											
	ID: DEC01 Decision type: Formula Predecessor: IPC Recurrences Exact: Min.: Max.: Rec.: Formula: if {Value activity Measured value activity MetaT002 / BO2 / IPC / HUM01} > {Measured value: max. intervention Measured value activity MetaT002 / BO2 / IPC / HUM01} then "DRY" else "NoDry" endif												
	Successors												
	<table><tr><th>ID</th><th>Condition</th><th>Description</th><th>Long description</th></tr><tr><td>DRY</td><td>DRY</td><td>Dry material</td><td></td></tr><tr><td>COMP</td><td>NoDry</td><td>No Drying Necessary</td><td></td></tr></table>	ID	Condition	Description	Long description	DRY	DRY	Dry material		COMP	NoDry	No Drying Necessary	
ID	Condition	Description	Long description										
DRY	DRY	Dry material											
COMP	NoDry	No Drying Necessary											
	Date / signature for basic function User identification: No identification User ident. group: ----- Verification signature: no												


1	(1) MetaT002/BO2/DRY Dry Material Active: yes Optional: no (Common BF)						
	Subtype: If the humidity is above 60%, dry material for max. 30 minutes						
1.0	REW01 - Dry Material Active: yes Subtype: None Done ---						
	10 Dry material in dryer D100 for max. 30 minutes						
	<table><tr><th>Set value</th><th>Actual value</th><th>Entry / Signature</th></tr><tr><td>Done</td><td></td><td></td></tr></table>	Set value	Actual value	Entry / Signature	Done		
Set value	Actual value	Entry / Signature					
Done							
	Date / signature for basic function User identification: No identification User ident. group: ----- Verification signature: no						

Active: yes **Reconciliation:** no
Start condition: Predecessors started

1	(1) MetaT002/BO2/COMP Compression Step Active: yes Optional: no (Common BF)	
	Subtype: Setup tablet press process parameter and proceed with compression	
1.0	ORD01 - Batch ID Active: yes Subtype: None	Batch ID ---
	10 The current running batch number	
	Set value Actual value Entry / Signature	
	Batch ID	
1.1	PRECOMP - Precompression Force Active: yes Subtype: None	1.0 (≥ 0.8 ; ≤ 1.2) kN ---
	10 Set precompression force	
	Set value Actual value Entry / Signature	
	1.0 (≥ 0.8 ; ≤ 1.2) kN	
1.2	MCF01 - Main Compression Force Active: yes Subtype: None	2.5 (≥ 2.0 ; ≤ 3.0) kN ---
	10 Set main ompression force	
	Set value Actual value Entry / Signature	
	2.5 (≥ 2.0 ; ≤ 3.0) kN	
1.3	TX01 - Program Number Active: yes Subtype: None	METAGRAM01 ---
	10 Set program number on compression equipment	
	Set value Actual value Entry / Signature	
	METAGRAM01	
1.4	Start - Compression start time Active: yes Subtype: None	---
	10 Please enter the timestamp of the start of compression	
	Set value Actual value Entry / Signature	
1.5	Color - Check Tablet Color Active: yes Subtype: None	Metagram Blue; Metagram White ---
	10 Select current tablet color from defined list	
	Set value Actual value Entry / Signature	
	Metagram Blue; Metagram White	
	Date / signature for basic function	
	User identification: No identification	
	User ident. group:	-----
	Verification signature: no	

1	(1) MetaT002/BO2/YS Yield & Stock Active: yes Optional: no (Common BF)
	Subtype: Determine yield and create stock Date / signature for basic function User identification: No identification User ident. group: ----- Verification signature: no

Active: yes Reconciliation: no
Start condition: Predecessors started



Page 13 of 14

Active: yes **Reconciliation:** no
Start condition: Predecessors started

L 2	(2) MetaT002/BO2/YS/DEC02 Create further HU Optional: NO (DECISION) Subtype: NONE												
	<table><tr><th>Actual value</th><th>Date</th><th>Signature</th></tr><tr><td></td><td></td><td></td></tr></table>	Actual value	Date	Signature									
Actual value	Date	Signature											
	If further Handling Units are needed, repeat step ID: DEC02 Decision type: User selection Predecessor: YD Recurrences Exact: Min.: Max.: Rec.: Successors <table><tr><th>ID</th><th>Condition</th><th>Description</th><th>Long description</th></tr><tr><td>MER02</td><td>ADDHU</td><td>ADDHU</td><td>Create further HU</td></tr><tr><td>SC</td><td>STOCK</td><td>STOCK</td><td>Create stock without further HU</td></tr></table> Date / signature for basic function User identification: No identification User ident. group: ----- Verification signature: no	ID	Condition	Description	Long description	MER02	ADDHU	ADDHU	Create further HU	SC	STOCK	STOCK	Create stock without further HU
ID	Condition	Description	Long description										
MER02	ADDHU	ADDHU	Create further HU										
SC	STOCK	STOCK	Create stock without further HU										

L 2	(2) MetaT002/BO2/YS/SC Stock Creation Active: yes Optional: no (Stock Creation) Subtype: None						
	<table><tr><th>Actual value</th><th>Date</th><th>Signature</th></tr><tr><td></td><td></td><td></td></tr></table>	Actual value	Date	Signature			
Actual value	Date	Signature					
	Create stock for all target material containers Target material: [100] 1 kg Meta0001 (Meta Aspirin) Attributes AQL: yes REM packing: no Stock creation type: Handling unit Execution mode: Manual						

L 2	(2) MetaT002/BO2/YS/YC Calculate Yield Active: yes Optional: no (Common BF) Subtype:						
	System calculates yield according to the actual produced quantity and batch size						
2.0	YC01 - Calculate Yield Active: yes Subtype: None <div>{MO Actual order size Master Batch Record MetaT002} / {MO Order size Master Batch Record MetaT002} * 100</div>						
	10 Execute formula for yield calculation						
	<table><tr><th>Set value</th><th>Actual value</th><th>Entry / Signature</th></tr><tr><td>{MO Actual order size Master Batch Record MetaT002} / {MO Order size Master Batch Record MetaT002} * 100</td><td></td><td></td></tr></table>	Set value	Actual value	Entry / Signature	{MO Actual order size Master Batch Record MetaT002} / {MO Order size Master Batch Record MetaT002} * 100		
Set value	Actual value	Entry / Signature					
{MO Actual order size Master Batch Record MetaT002} / {MO Order size Master Batch Record MetaT002} * 100							
2.1	END - Production Finished Active: yes Subtype: None						
	10 Please confirm the end date of the current production						
	<table><tr><th>Set value</th><th>Actual value</th><th>Entry / Signature</th></tr><tr><td></td><td></td><td></td></tr></table>	Set value	Actual value	Entry / Signature			
Set value	Actual value	Entry / Signature					
	Date / signature for basic function User identification: No identification User ident. group: ----- Verification signature: no						

1	(1) MetaT002/BO2/EQDAS1 EQM Disassembly Active: yes Optional: no (EQM Disassembly) Subtype: None						
	<table><tr><th>Actual value</th><th>Date</th><th>Signature</th></tr><tr><td></td><td></td><td></td></tr></table>	Actual value	Date	Signature			
Actual value	Date	Signature					
	Uninstall Punch Set from the tablet press Attributes Identification type: Scan Allowed status formula: Implicit release: no Multi-step EQM disass.: no Reuse equipment: no Equipment item: 1: Workroom MetaTab Meta Tableting						